# Characterization of human milk exosomes by Infrared Spectroscopy and LC-HRMS

Victoria Ramos-Garcia 1\*, Isabel Ten-Doménech1, Alba Moreno-Giménez1, Álvaro Solaz-García1, Inmaculada Lara-Cantón1, Amparo Ramón-Beltrán2, María Gormaz<sup>3</sup>, María Círia<sup>4</sup>, Pilar Sepúlveda<sup>4</sup>, David Pérez-Guaita<sup>5</sup>, Bernhard Lendl<sup>6</sup>, Guillermo Quintás<sup>7,8</sup>, Julia Kuligowski<sup>1</sup>

<sup>1</sup>Neonatal Research Unit, Health Research Institute Hospital La Fe, Avda Fernando Abril Martorell 106, 46026 Valencia, Spain

<sup>2</sup>Human Milk Bank, University & Polytechnic Hospital La Fe, Avda Fernando Abril Martorell 106, 46026 Valencia, Spain

<sup>3</sup>Division of Neonatology, University & Polytechnic Hospital La Fe, Avda Fernando Abril Martorell 106, 46026 Valencia, Spain 4Regenerative Medicine and Heart Transplantation Unit. Health Research Institute Hospital La Fe. Avda Fernando Abril Martorell 106. 46026 Valencia. Spain

<sup>5</sup>Department of Analytical Chemistry, University of Valencia, 50 Dr. Moliner Street, research building, 46100 Burjassot, Valencia, Spain

<sup>6</sup>Institute of Chemical Technologies and Analytics, Technische Universität Wien, Getreidemarkt 9/164. A 1060 Vienna, Austria

<sup>7</sup>Health and Biomedicine, Leitat Technological Center, Carrer de la Innovació, 2, 08225 Terrassa, Spain

<sup>8</sup>Analytical Unit, Health Research Institute La Fe, Ayda Fernando Abril Martorell 106, 46026 Valencia, Spain

\*victoria\_ramos@iislafe.es

# **OBJECTIVES**

La Fe

Instituto de Investigación

Sanitaria La Fe

Universitari

i Politècnic

Exosomes are nanosized vesicles containing specific cargos of DNA, RNA, proteins, metabolites, and intracellular and membrane lipids. Lipid analysis of exosomes is a relatively unexplored field of research, so the main objective of this study was to:

• Determine the feasibility of Attenuated Total Reflectance - Fourier Transform Infrared (ATR-FTIR) spectroscopy for the routine control of isolated exosomes and the rapid characterization of their lipid profiles.


# **EXOSOME ANALYSIS**



- ✓ Bicinchoninic acid (BCA) assay to determine the amount of retrieved exosomes
- ExoView platform (NanoView Biosciences, MA, USA) to determine the exosome size distribution and quantification of vesicles

<sup>1</sup>Ten-Doménech, et al., Comparing Targeted vs. Untargeted MS2 Data-Dependent Acquisition for Peak Annotation in LC-MS Metabolomics. Metabolites 2020, 10, 126.

## STUDY DESIGN



Multivariate analysis METHIN HIMOD

LipidBlast

## ATR-FTIR

- √ 2 µL dried HM exosomes were used for IR spectra acquisition
- Range: 4000 to 400 cm<sup>-1</sup>
- Scans: 32
- Resolution: 4 cm<sup>-1</sup>




metabolites extract


# **UPLC-QqTOF-MS**

- ✓ Untargeted HM and HM exosomes analysis by ultra-high performance liquid chromatography quadrupole-time-of-flight mass spectrometry (UPLC-qTOF-MS) operating in positive and negative ionization mode (ESI+ and ESI-)
- Column: Acquity BEH C18 (100 mm x 2.1 mm, 1.7μm)
- Mobile phase: (A) 5:1:4 IPA:CH<sub>3</sub>OH:H<sub>2</sub>O 5mM CH<sub>3</sub>COONH<sub>4</sub> 0.1%HCOOH; (B) 99:1 IPA:H<sub>2</sub>O 5mM CH<sub>3</sub>COONH<sub>4</sub> 0.1%HCOOH
- Automated MSMS-based annotation of metabolites<sup>1</sup>

# **QUALITY CONTROL**




# **QUALITATIVE ANALYSIS**



Statistical heterospectroscopy was employed to assess the correlation among IR bands and LC-MS intensities of the annotated features showing significant associations

IR spectra (top) of dry residues of exosomes and the first (middle) and second derivative (bottom). Black line: PBS blank




A) Principal component analysis was used to reveal samples with similar lipidomic profiles, as they were clustered mainly along PC1 (p < 0.05). B) and C) Sum of intensities of lipids annotated in the three selected clusters for LC-MS and second derivative IR spectra, respectively

# **RESULTS**

# **UPLC-QqTOF-MS** 377 features annotated in the data set after data pre-processing and clean-up

# QUANTITATIVE ANALYSIS



Quantitative analysis using A) the area of the Amide I band and the results of the BCA assay;  ${\bf B}{\bf )}$  the same band as in A) of an albumin external calibration and the results of the BCA assay; C) and D) the sum of LC-MS intensities of annotated lipids as surrogate of the total lipid content and the area of the CHs and C=O streching bands, respectively

- A correlation between the lipidomic profile of HM exosomes and their ATR-FTIR spectra has been obtained
- The latter technique can be used for a rapid evaluation of the composition of HM exosomes (both, qualitative and quantitative)
- A new tool for a direct and fast quality control of the exosome isolation procedure

### ACKNOWLEDGEMENTS